

## PCR CLIENT SERVICES MANUAL

www.streamlinesci.com

## WELCOME LETTER

Dear Customer,

Thank you for trusting Streamline Scientific for your molecular, microbiology, and toxicology testing. As a physician-founded company, we are passionate about the services we provide and know that the appropriate diagnostic information informs your decision making, and improves patient care.

Since our start in 2011, we have focused on simplifying the delivery of molecular technologies for both medical practices and patients. We do this by arming clinicians with the insights needed to prescribe the most appropriate treatment - within 24 hours.

Our sister organization, Molecular Designs, is a leading developer and manufacturer of PCR assays. In March of 2020, upon the news of the novel coronavirus, Molecular Designs was one of the first of 5 labs in the nation to submit an EUA and be granted FDA authorization.

This spirit of innovation continues today, developing novel tests to meet market needs and exploring new ways to better serve our many Alabama-based customers. Several of the elements that define the "Streamline Difference" include:

- Industry-leading technology and customer portal
- EMR integrations

- Unmatched customer support and dedicated team
- Door-to-door courier service
- An improved patient experience with multiple tests from a single swab

You are in exceptional hands with our Business Development and Client Relations teams, and the enclosed documentation will serve as a useful resource.

Should you have any questions or input on how we can better serve you, please don't hesitate to contact me directly.

Sincerely yours,



SHAWN HOOD

President shood@StreamlineSci.com



## DIRECTORY

## BIRMINGHAM

2868 Acton Road Suite 207 Birmingham, AL 35243

CLIA - 01D2074949 COLA - 25348



## AUSTIN

13413 Galleria Circle Suite Q-140 Bee Cave, TX 78738

CLIA – 45D2268534 COLA – 32244

## CLIENT SERVICES CONTACT INFORMATION

Monday - Sunday: 8AM - 5PM

Email: clientservices@streamlinesci.com

Local: 205.558.0095

Toll Free: 855.319.4459

Fax: 877.796.6185

## BILLING SERVICES CONTACT INFORMATION

Email:	patientbilling@streamlinesci.com
--------	----------------------------------

Local: 844.364.9969

## LAB OPERATIONS

7 days a week

(Any holiday closures will be announced in advance and posted on streamlinesci.com)



## SAMPLE COLLECTION GUIDELINES

### **COLLECTION CONTAINERS**

- Synthetic Flocked Swab with Liquid Amies Transport buffer
- Biotech Saline Swab (Nasopharyngeal/ Oropharyngeal) – COVID and Respiratory Panels ONLY
- Urine Vacutainers
- Sterile Toxicology Cup with temperature strip (labeled as: Clicktainer Vial, Temperature Strip Label, Assembled) – Toxicology Samples Only
- Sterile Dry Tubes (labeled as Transport Tube, 5ml, Capped, Self-Standing,

Sterile, case 50) – Fungal Infection Panel ONLY

### SPECIMEN TYPE CONSIDERATIONS

#### PCR Testing:

- Swabs Ideal samples are taken from the source of the potential infection.
   Ex. Nasopharyngeal, oropharyngeal, throat, superficial/subdermal wounds, vaginal tract, etc.
- Urine If sending in urine specimens, clean-catch urine is required for UTI.
   Clean-catch urine is not recommended for STI or Candida testing.

#### Antimicrobial Susceptibility Testing:

- Swabs Can be used from most sources. For pure nasal swabs, only culturing performed (Antimicrobial Susceptibility Testing is not recommended for nasal swabs due to the lack of clinical significance).
- Urine (bacteriostatic preservative) No special requirements.

### **REJECTION CRITERIA**

#### General criteria (applicable to all tests and sample types):

- Specimen is unlabeled or improperly labeled. Properly labeled includes two patient identifiers matching submitted requsition form.
- Requisition form is missing or incomplete.
- Specimen is improperly sealed, resulting in leakage and possible contamination during transport.
- Incorrect specimen type for selected panel.
- Insufficient sample volume for testing.

#### Swab Specimens:

- Swab specimens received with no original swab present in tube.
- Dry swab with no buffer present in tube.

#### Urine Specimens for UTI Testing:

• Urine specimens submitted in Aptima buffer or other container with no **bacteriostatic** preservative.

#### PCR Testing (any panel):

Specimens received after 5 days post-collection.

#### Culturing and Antimicrobial Susceptibility Testing:

- Swab specimens received after 2 days post-collection.
- Swab specimens received in a non-liquid Amies buffer (such as viral transport medium).
- Nasal swab specimens (exceptions see specimen type considerations).
- Urine specimens received after 3 days post-collection.
- STI or Vaginitis specimens Antimicrobial Susceptibility Testing is not available.

#### **GENERAL DISCLAIMERS**

**Applies to all patient specimens submitted to Streamline Scientific:** Proper patient identification and accurate specimen labeling is required per CLIA (Clinical Laboratory Improvement Amendment). Streamline Scientific has determined the analytical performance characteristics of the molecular tests. They have not been cleared or approved by the U.S. Food and Drug Administration. It is preferable that all specimen types are received within 2 days of collection. They will need to be received within the specimen stability range stated by the manufacturer of the collection device measured from the time of collection to time of receipt for laboratory testing, to prevent a disclaimer on the test report.



Alabama - CLIA: 01D2074949 Texas - CLIA: 45D2268534 **VERSION: 1.1.23** 🔊 streamline 2868 Acton Road, Suite 207, Birmingham, AL 35243 13413 Galleria Circle, Suite Q-140, Bee Cave, TX 78738 p: 855.319.4459 | f: 877.796.6185 | StreamlineSci.com p: 855.319.4459 | f: 877.796.6185 | StreamlineSci.com Male Last Name / First Name / Client Name / Account Asian Address / Address / APT# Black **STEP 1 STEP 2** City / State / Zi Native American Other **Ordering Provider** Subscriber ID Non-Hispanic Bill to: 🗌 Insurance 🗌 Facility Group # necessary to rapidly determine d timely data available **STEP 3 STEP 4** use (50% according provider with my for all co-pays an Verbal Order **Provider Signature** Patient Signature Standing Order COVID-19 Only UTI w/ ABX Resistance Wound/Derm w/ ABX Fungal Infection COVID Respiratory Lite (includes all the pathogens in the pane Haemophilus influenzae Treundii ) Citrobacter koseri ) Enterobacter cloacae ) Enterococcus spp. ) Escherichia coli **STEP 5** N76.0 Acute vaginitis N89.8 Other specified noninflammatory disorders of vagina R36.9 Urethral discharge unspecified Morganella morganii Proteus mirabilis Pseudomonas aeruginosa Z30.9 Encounter for contraceptive management Trichophyton zoophilic spp Microsporum canis Bacterial Add On ABX Resistance Marker Methicillin/Oxacillin (mecA) ABX Resistance Markers ABX Resistance Markers SPECIMEN SOURCE UTI Plus Z22.39 Carrier of other specified bacterial diseases N77.1 Vaginitis, vulvitis, & vulvovaginitis B37.3 Candidiasis of vulva & vagina Candida dubliniensis Z22.322 Carrier or suspected carrier of MRSA ) Staphylococcus epidermidis ) ABX Resistance Marker Candida glabrata Candida krusei Z30.9 Encounter for contraceptive management ICD 10 CODES R05.9 Cough, unspecified N76.0 Acute N89.8 Other specified noninflammatory ICD 10 CODES Culture ID w/ Reflexive Antimicrobial Susceptibility Testing (AST is not available for STI or Vaginitis) SPECIMEN SOURCE Rectal Swab Other 228.310 Partially vascinated for COVID-228.310 Unvascinated for COVID-19 228.39 Other underimmunization status SPECIMEN SOURCE Nasal Swab Ear Swab Nasopharyngeal Swab SPECIMEN SOURCE **STEP 6** Clean catch urine 6 Please indicate if your patient has taken antibiotics in tr No

FOR MOST UP TO DATE REQUISITION FORM, PLEASE VISIT WWW.STREAMLINESCI.COM/FORMS-SUMMARY-LINKS



## REQUISITION FORM INSTRUCTIONS

### **STEP 1**

#### PATIENT INFORMATION

Fill in the patient's information including their name, gender, address, race, phone number, email, date of birth, social security number, insurance, ethnicity and billing information.

## STEP 2

#### **PROVIDER INFORMATION**

Fill in the provider's information, including client name or account number, address, phone number, ordering physician, collection date, specimen collector, collection time, and state where sample was collected.

### **STEP 3**

#### MEDICAL NECESSITY

Have the physician read, sign, and choose either a verbal or standing order for the medical necessity.

## STEP 4

#### CONSENT FOR TESTING

Have the patient read and sign the medical necessity.



#### PANEL LIST

Check appropriate panels that address your patients' needs. Tests can be ordered individually. Select the box(es) that describes the specimen source.

Make sure to check any ICD-10 codes found at the end of each panel.



Patient Information

Address /

Phone

DOB

Insurance

Group #

Provider Information

Address / APT#

City / State / Zip

Ordering Provider

Specimen Collecte

3 Medical Necessity

Provider Signature

Provider Information

Patient Signature

Phone

Client Name / Account

Last Name / First Name /

City / State / Zip / County

Email

SSN

Subscriber ID Bill to: Insurance Facility

Fax #

State Co

As part of my antibiotic stewardship policy, I find it medically necessary to rapidly determine and differentiate a viral and/or bacterial infection in order to treat with or without appropriate antibiotics. Having the most accurate and timely data avait to me directly guides my treatment and patient management. Empiric treatment a

The information I have provided on this form is accurate. I authorize Streamline Scientific to release the results of this test to my treating physician of for services I receive I am aware that Streamline Scientific may be an out of exit provider with my insurer. I am aware that I am responsible for all copys and deductibles not covered by insurance or other payers.

Date

management leads to inappropriate and unnecessary to the CDC) and delayed diagnosis which can lead to □Male □Female

Race: Asian Black

Black Caucasian Hispanic Native American Other N/A

Ethnicity: Hispanic Non-Hispanic

Collection Date

Collection Time

Verbal Order

### STEP 6

#### HAS PATIENT TAKEN ANTIBIOTICS? (3 Please indicate if your patient has taken antibiotics in the past 72 hours:

Check "yes" if your patient has taken antibiotics within the past 72 hours or "no" if they have not.

🗌 Yes 🗌 No



## SPECIMEN PACKAGING, TRANSPORT, AND RESULT DELIVERY



## SPECIMEN PACKAGING AND TRANSPORT

- 1. Ensure sample lids are secure to avoid leakage or contamination
- 2. Place specimen into a provided biohazard bag
- 3. Insert requisition form, including all patient insurance information and demographics, into the specimen biohazard bag's side pocket
- 4. Place the biohazard bag containing the patient specimen into the provided lock box.
- 5. Enter a pickup request ticket through the online portal or contact Streamline Scientific customer service team to schedule a pickup (unless a nightly pickup schedule has been established).





## RESPIRATORY PANEL TARGET CATALOG

## VIRAL PATHOGENS

- Adenovirus
- Bocavirus
- COVID-19 (SARS-CoV-2)
- Coronavirus 229E
- Coronavirus HKU1
- Coronavirus NL63
- Coronavirus OC43
- EBV (mononucleosis)
- Enterovirus
- Human Metapneumovirus A
- Human Metapneumovirus B
- Influenza A
- Influenza B
- Parainfluenza Virus type 1
- Parainfluenza Virus type 2
- Parainfluenza Virus type 3
- Parainfluenza Virus type 4
- Respiratory Syncytial Virus
- Rhinovirus

## BACTERIAL PATHOGENS

- Acinetobacter baumannii\*
- Bordetella pertussis
- Chlamydophila pneumoniae
- Enterobacter cloacae\*
- Haemophilus influenzae\*
- Klebsiella aerogenes\*
- Klebsiella pneumoniae\*
- Legionella pneumonophila
- Moraxella catarrhalis\*
- Mycoplasma pneumoniae
- Proteus mirabilis\*
- Pseudomonas aeruginosa\*
- Staphylococcus aureus\*
- Staphylococcus epidermidis
- Streptococcus pneumoniae\*
- Streptococcus pyogenes (Group A)\*\*

## ABX RESISTANCE MARKER

• Methicillin/Oxacillin (mecA)

\*Limitation: Test does not differentiate between a patient with acute infection or an asymptomatic carrier.

\*\*S. pyogenes detected from a throat swab is diagnostic of pharyngitis; S. pyogenes detected from a nasopharyngeal swab could indicate an asymptomatic carrier.

For panel offerings, please visit streamlinesci.com



## RESPIRATORY INFECTION SAMPLE COLLECTION

Follow the instructions below depending on the specimen source:









### NASOPHARYNGEAL SWAB:

- 1. Tilt the patient's head back 70 degrees.
- 2. Insert the swab into the nostril (the swab should reach depth equal to the distance from nostril to outer opening of the ear). Leave the swab in place for several seconds to absorb secretions.
- 3. Slowly remove the swab while rotating it.

### MID-TURBINATE (PEDIATRIC) SWAB:

- 1. Gently insert the swab into the nostril.
- 2. Using a gentle rotation, push the swab until a slight resistance is met at the level of turbinates.
- 3. Rotate the swab several times against the nasal wall.

### **OROPHARYNGEAL SWAB:**

- 1. Insert the swab into the posterior pharynx and tonsillar areas.
- 2. Rub the swab over both tonsillar pillars and posterior oropharyngeal and avoid touching the tongue, teeth, and gums.

### **SPUTUM:**

- 1. Educate the patient about the difference between sputum and oral secretions
- 2. Have the patient rinse the mouth with water and then expectorate deep cough sputum directly into a sterile screw-cap collection cup or sterile dry container.
- 3. Transfer this sputum to a Copan E-swab tube with a liquid buffer for transport







2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION PATIENT INFORMATION SPECIMEN INFORMATION** Name: Streamline Pathology PATIENT, TEST Name: Lab Accession Number: S-08232210321 Laboratory 6/30/1980 DOB: 08/22/22 12:00:00 AM Sex: M Date Collected: Address: 2868 Acton Road Phone: (205)332-3160 Date Accessioned: 08/23/2022 Suite 207 BIRMINGHAM, AL 35243 2868 ACTON ROAD Address: Date Reported: 08/23/2022 Provider: Test, Doctor MD VESTAVIA, AL35243 Faxed to: 18777966185

#### Controls

Patient Extraction Control <sup>1</sup> PASS Endogenous Positive Control <sup>1</sup> PASS Pathogen Positive Control <sup>2</sup> PASS Pathogen Negative Control <sup>3</sup> PASS

(1) Endogenous control confirms sample collection, DNA/RNA extraction, and assay enzyme activity (2) Positive control is synthetic inactive pathogen (3) Negative Control contains primers, probe, and enzymes with no DNA/RNA template (4) A "Detected" result indicates the presence of a pathogen (99.99% confidence) above the assay cutoff.

**Test Performed** 

**Test Results** 

Comments

COVID-19		Collection Type: Nasopharyngeal Swab		
COVID-19	Not Detected	[08/23/22] COVID-19 assay reviewed and approved under FDA Emergency Use Authorization #200522.		
nCoV_N1 (CT Value)	0	[08/23/22] CT value indicates the number of amplification cycles of real-time PCR needed to detect COVID-19, a specific SARS-CoV-2 gene sequence. CT values are inversely proportional to the amount of target nucleic acid in the sample (i.e. the lower the CT level, the greater the amount of target nucleic acid in the sample). The reference range for Assurance Scientific Laboratories RT-PCR assays is <40 cycles of amplification.		

Test Performed	Lab Result (4)	DNA Сору	Commonte
Test Performed	(Qualitative Low/Medium/High)	Number	comments

SUMMARY COVID Respiratory	/ Plus	Collection Type: Nasopharyngeal Swab	
Acinetobacter baumannii DETECTED - HIGH		1.00E+06	[08/23/22] Limitation: Test does not differentiate between patient with acute infection or an asymptomatic carrier.
mecA (Methicillin/Oxacillin resistance)	DETECTED	DETECTED	
Rhinovirus (types A & B)	DETECTED - MEDIUM	1.00E+03	[08/23/22] Assay is developed to detect all strains of this pathogen
			[08/23/22] May cross-react with Enterovirus
Streptococcus pneumoniae	DETECTED - LOW	1.00E+01	[08/23/22] Limitation: Test does not differentiate between patient with acute infection or an asymptomatic carrier.

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Page: 1

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE RESPIRATORY PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com

streamline scientific



2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### CLINIC INFORMATION

Name:	Streamline Pathology Laboratory	Na
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243	Pł
Provider:	Test, Doctor MD	A

PATIENT INFORMATION					
Name:	PATIENT, TEST				
DOB:	6/30/1980	Sex:	М		
Phone:	(205)332-3160				
Address:	2868 ACTON ROAD VESTAVIA, AL35243				

#### SPECIMEN INFORMATION

Lab Accession Number:	S-08232210321
Date Collected:	08/22/22 12:00:00 AM
Date Accessioned:	08/23/2022
Date Reported:	08/23/2022
Faxed to:	18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments	
COVID Respiratory Plus			Collection Type: Nasopharyngeal Swab	
Acinetobacter baumannii	DETECTED - HIGH	1.00E+06	[08/23/22] Limitation: Test does not differentiate between patient with acute infection or an asymptomatic carrier.	
Adenovirus	Not Detected	Not Detected	[08/23/22] Assay is developed to detect all strains of this pathogen	
Bocavirus	Not Detected	Not Detected	[08/23/22] Assay is developed to detect all strains of this pathogen	
Bordetella pertussis	Not Detected	Not Detected	[08/23/22] Shares homology with B.holmesii.	
Chlamydophila pneumoniae	Not Detected	Not Detected		
Coronavirus 229E	Not Detected	Not Detected		
Coronavirus HKU1	Not Detected	Not Detected		
Coronavirus NL63	Not Detected	Not Detected		
Coronavirus OC43	Not Detected	Not Detected		
EBV (Mononucleosis)	Not Detected	Not Detected		
Enterobacter cloacae	Not Detected	Not Detected		
Enterovirus	Not Detected	Not Detected	[08/23/22] Enterovirus includes Coxsackievirus types A9, A10, A16, B5 and Echovirus serotypes	
			[08/23/22] May cross-react with Rhinovirus (types A and B)	
Haemophilus influenzae	Not Detected	Not Detected		
HMPV A (Human Metapneumovirus)	Not Detected	Not Detected		
HMPV B (Human Metapneumovirus)	Not Detected	Not Detected		
Influenza A	Not Detected	Not Detected		
Influenza B	Not Detected	Not Detected		
Klebsiella aerogenes	Not Detected	Not Detected		
Klebsiella pneumoniae	Not Detected	Not Detected		
Legionella pneumophila	Not Detected	Not Detected		
mecA (Methicillin/Oxacillin resistance)	DETECTED	DETECTED		
Moraxella catarrhalis	Not Detected	Not Detected		
Mycoplasma pneumoniae	Not Detected	Not Detected		
Parainfluenza 1	Not Detected	Not Detected		
Parainfluenza 2	Not Detected	Not Detected		
Parainfluenza 3	Not Detected	Not Detected		
Parainfluenza 4	Not Detected	Not Detected		
Proteus mirabilis	Not Detected	Not Detected		
Pseudomonas aeruginosa	Not Detected	Not Detected		

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE RESPIRATORY PANEL SAMPLE LAB REPORTS, CLICK HERE



Page: 2



2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Name:	Streamline Pathology Laboratory
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243
Provider:	Test, Doctor MD

PATIENT INFORMATION					
Name:	PATIENT, TEST				
DOB:	6/30/1980	Sex:	М		
Phone:	(205)332-3160				
Address:	2868 ACTON ROAD VESTAVIA, AL35243				

# SPECIMEN INFORMATIONLab Accession Number:S-08232210321Date Collected:08/22/22 12:00:00 AMDate Accessioned:08/23/2022Date Reported:08/23/2022Faxed to:18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Rhinovirus (types A & B)	DETECTED - MEDIUM	1.00E+03	[08/23/22] Assay is developed to detect all strains of this pathogen
			[08/23/22] May cross-react with Enterovirus
RSV A/B (Respiratory Syncytial Virus)	Not Detected	Not Detected	
Staphylococcus aureus	Not Detected	Not Detected	
Staphylococcus epidermidis	Not Detected	Not Detected	
Streptococcus pneumoniae	DETECTED - LOW	1.00E+01	[08/23/22] Limitation: Test does not differentiate between patient with acute infection or an asymptomatic carrier.
Streptococcus pyogenes (Group A)	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 3

streamline

scientific

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE RESPIRATORY PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com

## UTI PANEL TARGETS

## PATHOGENS

- Acinetobacter baumannii
- Bacteroides fragilis
- Candida albicans
- Candida dubliniensis
- Candida glabrata
- Candida krusei
- Candida parapsilosis
- Candida tropicalis
- Citrobacter braakii/freundii
- Citrobacter koseri
- Enterobacter cloacae
- Enterococcus spp.
- Escherichia coli
- Klebsiella aerogenes
- K. oxotyca/michiganensis

- Klebsiella pneumoniae
- Morganella morganii
- Mycoplasma genitalium
- Mycoplasma hominis
- Prevotella bivia
- Proteus mirabilis
- Pseudomonas aeruginosa
- Serratia marcescens
- Staphylococcus aureus
- Staphylococcus epidermidis
- Staphylococcus saprophyticus
- Streptococcus agalactiae (Group B)
- Streptococcus pyogenes (Group A)
- Ureaplasma urealyticum

## **ABX RESISTANCE MARKER**

- β-lactamase (blaKPC)
- β-lactamase (CTX-M-Group 1)
- metallo-β-lactamase (blaNDM)
- Fluoroquinolones

- Methicillin/Oxacillin (mecA)
- Sulfonamides
- Trimethoprim

For panel offerings, please visit streamlinesci.com





## UTI SAMPLE COLLECTION

Collect the patient sample using one of the procedures below:

### **URINE SAMPLE:**

If submitting a urine sample: Collect a urine sample from the patient in a urine collection cup with lid-integrated transfer device using the following gender-specific instructions:

### **Female Urine Sample**

- 1. Patient should use a packaged, moist towel to clean the vulva and perianal area starting from front to back. Repeat with a second moist towel.
- 2. Patient should then spread their labia with one hand and start urinating into the toilet. With the other hand, they should put the urine container under the genital area to catch the stream of urine without touching any skin.

### **Male Urine Sample**

- 1. Patient should retract the foreskin from the penis if necessary and use the packaged towel to clean the penis from the tip to the base. Repeat with a second moist towel.
- 2. Patient should retract the foreskin if necessary with one hand and start urinating into the toilet. Then, position the urine conainer with the other hand to catch the stream without touching any skin.

Replace and tighten the collection cup lid, and transfer the urine sample into a vacutainer tube with preservative:

- 1. Remove the sticker from the lid to access the integrated transfer device.
- 2. Insert the vacutainer tube vertically into the transfer device, puncturing the rubber septum and allowing the tube to fill with urine.
- 3. Remove the vacutainer tube and discard the collection cup. Mix urine and preservative by gently inverting the tube 8-10 times.











2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

CLINIC INFORMATION PA		PATIENT	PATIENT INFORMATION		SPECIMEN INFORMAT	ΓΙΟΝ	
Name:       Streamline Pathology Laboratory       Na Du Du Du Du Suite 207 BIRMINGHAM, AL 35243         Provider:       Test, Doctor MD		Name: DOB: Phone: Address:	PATIENT, TEST 6/30/1980 Sex: M (205)332-3160 2868 ACTON ROAD VESTAVIA, AL35243 (1) Endogenous control confirms sampl (2) Positive control is synthetic ina		Lab Accession Number: Date Collected: Date Accessioned: Date Reported: Faxed to:	S-08232210334 08/22/22 12:00:00 AM 08/23/2022 08/23/2022 18777966185	
Endogenous Positive Control <sup>1</sup> PASS Pathogen Positive Control <sup>2</sup> PASS Pathogen Negative Control <sup>3</sup> PASS		1   PASS     2   PASS     3   PASS	<pre>(3) Negative (4) A "Detec cutoff.</pre>	e Control con	tains primers	probe, and enzymes with no E presence of a pathogen (99.99	NA/RNA template % confidence) above the assay
Test Pe	rformed	Lab Result (Qualitative Low/M	t (4) edium/High)	DNA Co Numb	opy Jer	Comme	ents
SUMMAR	RY UTI Pathogens Plus				Co	ollection Type: Urine clear	n catch
Acineto	bacter baumannii	DETECTED	- HIGH	1.00E+	05		
Morgan	ella morganii	DETECTED -	MEDIUM	1.00E+	.03		
Streptor	coccus agalactiae (GBS)	DETECTED	- LOW	1.00E+	02		

Class (Gene Name)	Lab Result (Qualitative)	Resistance Gene Targets Identified	Associated Resistances (Antibiotics to Avoid)
SUMMARY UTI Antibiotic Resist	ance Markers		Collection Type: Urine clean catch
Class A Beta-lactamase (CTX-M- Group 1)	DETECTED	blaCTX-M-1,3,10,12,15,22,23,28; blaFEC-1	[08/23/22]Cephalosporins, Penicillins, Aztreonam

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE UTI PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com

streamline scientific

Page: 1



М

2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Name:	Streamline Pathology Laboratory
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243
Provider:	Test, Doctor MD

PATIENT	INFORMATION	
Name:	PATIENT, TEST	
DOB:	6/30/1980	Sex
Phone:	(205)332-3160	
Address:	2868 ACTON ROAD VESTAVIA, AL35243	

#### SPECIMEN INFORMATION

Lab Accession Number:	S-08232210334
Date Collected:	08/22/22 12:00:00 AM
Date Accessioned:	08/23/2022
Date Reported:	08/23/2022
Faxed to:	18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
UTI Pathogens Plus			Collection Type: Urine clean catch
Acinetobacter baumannii	DETECTED - HIGH	1.00E+05	
Bacteroides fragilis	Not Detected	Not Detected	
Candida albicans	Not Detected	Not Detected	
Candida dubliniensis	Not Detected	Not Detected	
Candida glabrata	Not Detected	Not Detected	
Candida krusei	Not Detected	Not Detected	
Candida parapsilosis	Not Detected	Not Detected	
Candida tropicalis	Not Detected	Not Detected	
Citrobacter freundii/braakii	Not Detected	Not Detected	
Citrobacter koseri	Not Detected	Not Detected	
Enterobacter cloacae	Not Detected	Not Detected	
Enterococcus spp.	Not Detected	Not Detected	
Escherichia coli	Not Detected	Not Detected	[08/23/22] Shares homology with Shigella spp.
Klebsiella aerogenes	Not Detected	Not Detected	
Klebsiella oxytoca/michiganesis	Not Detected	Not Detected	
Klebsiella pneumoniae	Not Detected	Not Detected	
Morganella morganii	DETECTED - MEDIUM	1.00E+03	
Mycoplasma genitalium	Not Detected	Not Detected	
Mycoplasma hominis	Not Detected	Not Detected	
Prevotella bivia	Not Detected	Not Detected	
Proteus mirabilis	Not Detected	Not Detected	
Pseudomonas aeruginosa	Not Detected	Not Detected	
Serratia marcescens	Not Detected	Not Detected	
Staphylococcus aureus	Not Detected	Not Detected	
Staphylococcus epidermidis	Not Detected	Not Detected	
Staphylococcus saprophyticus	Not Detected	Not Detected	
Streptococcus agalactiae (GBS)	DETECTED - LOW	1.00E+02	
Streptococcus pyogenes (Group A)	Not Detected	Not Detected	
Ureaplasma urealvticum	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Page: 2

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE UTI PANEL SAMPLE LAB REPORTS, CLICK HERE





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Trimethoprim

vanA (Vancomycin) vanB (Vancomycin)

Name:	Streamline Pathology Laboratory	Na	
Address:	2868 Acton Road	DO	
	Suite 207	Ph	
	BIRMINGHAM, AL 35243	Ad	
Provider:	Test, Doctor MD	1.0	

PATIENT INFORMATION				
Name:	PATIENT, TEST			
DOB:	6/30/1980	Sex:	Μ	
Phone:	(205)332-3160			
Address:	2868 ACTON ROAD VESTAVIA, AL35243			

#### **SPECIMEN INFORMATION**

Lab Accession Number:	S-08232210334
Date Collected:	08/22/22 12:00:00 AM
Date Accessioned:	08/23/2022
Date Reported:	08/23/2022
Faxed to:	18777966185

Class (Gene Name) Lab Result (Qualitative)		Resistance Gene Targets Identified	Associated Resistances (Antibiotics to Avoid)
UTI Antibiotic Resistance Marke	rs		Collection Type: Urine clean catch
Class A Beta-lactamase (blaKPC)	Not Detected	KPC-2-8,10,11,13-22,24-33	[08/23/22] Carbapenems, Cephalosporins, Penicillins, Beta-lactamase inhibitors, Aztreonam
Class A Beta-lactamase (CTX-M- Group 1)	DETECTED	blaCTX-M-1,3,10,12,15,22,23,28; blaFEC-1	[08/23/22] Cephalosporins, Penicillins, Aztreonam
Class B metallo Beta-lactamase (blaNDM)	Not Detected	NDM (1-21)	[08/23/22] Carbapenems, Cephalosporins, Penicillins, Beta-lactamase inhibitors
Fluoroquinolones	Not Detected	qnrS 1-5,7-9; qnrB Group 1; qnrB Group 5	[08/23/22] Ciprofloxacin, Gemifloxacin, Levofloxacin, Moxifloxacin, Norfloxacin, Ofloxacin
mecA (Methicillin/Oxacillin resistance)	Not Detected	mecA	[08/23/22] Oxacillin
Sulfonamides	Not Detected	sul1; sul2; sul3	[08/23/22] Sulfadiazine, Sulfamethizole, Sulfamethoxazole, Sulfasalazine, Sulfisoxazole
Trimethoprim	Not Detected	dfrA1; dfrA5; dfrA11; dfrA17	[08/23/22] Primsol

vanA

vanB

Not Detected

Not Detected

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

[08/23/22] Vancomycin

[08/23/22] Vancomycin

#### Lab Director: Ty Thomas, MD

Page: 3

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

FOR MORE UTI PANEL SAMPLE LAB REPORTS, CLICK HERE



## WOUND PANEL

## PATHOGENS

- Acinetobacter baumannii
- Bacteroides fragilis
- Citrobacter braakii/freundii
- Citrobacter koseri
- Enterobacter cloacae
- Enterococcus spp.
- Escherichia coli
- Klebsiella aerogenes
- K. oxytoca/michiganensis

- Klabsiella pneumoniae
- Morganella morganii
- Proteus mirabilis
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Staphylococcus epidermidis
- Staphylococcus saprophyticus
- Streptococcus pyogenes (Group A)
- Varicella Zoster (Shingles)

## **ABX RESISTANCE MARKER**

- β-lactamase (blaKPC)
- β-lactamase (CTX-M-Group 1)
- metallo-β-lactamase (blaNDM)
- Fluoroquinolones

- Methicillin/Oxacillin (mecA)
- Sulfonamides
- Trimethoprim

For panel offerings, please visit streamlinesci.com



## WOUND/INFECTION SAMPLE COLLECTION

Collect the patient sample using one of the procedures below:

### **ASPIRATION**

- The surface of the wound/abscess should be carefully cleansed and debrided using sterile gauze and saline before attempting to aspirate the specimen.
- 2. Aspirate the specimen and place 0.5 to 1.0 mL of the aspirate directly into a Copan eSwab® transport tube.
- 3. Without contaminating the included swab, place the swab into the Copan eSwab® transport tube all the way to the bottom. Break the swab at the scored breakpoint indication line and leave the bottom portion inside the transport tube, partially submerged in buffer solution. Screw the top onto the tube tightly to secure the specimen.

### **SWAB**

- 1. Cleanse and debride the wound with sterile gauze and saline.
- 2. Exudate and brushings of the wound base (including advancing margins) should be collected using the swab.
- 3. Without contaminating the swab, place the swab into the transport tube all the way to the bottom.
- 4. Break the swab at the scored breakpoint indication line and leave the bottom portion inside the Copan eSwab® transport tube, partially submerged in buffer solution. Screw the top onto the tube tightly to secure the specimen.









2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

CLINIC INFORMATION PATIENT INFORMAT			ΓΙΟΝ		SPECIMEN INFORMAT	ΓΙΟΝ	
Name: Address: Provider:	Streamline Pathology Laboratory 2868 Acton Road Suite 207 BIRMINGHAM, AL 352 Test, Doctor MD	Name: DOB: Phone: 43 Address:	PATIENT, 1 6/30/1980 (205)332-3 <sup>-1</sup> 2868 ACTO VESTAVIA,	EST Se 160 N ROAD AL35243	ex: M	Lab Accession Number: Date Collected: Date Accessioned: Date Reported: Faxed to:	S-08232210336 08/22/22 12:00:00 AM 08/23/2022 08/23/2022 18777966185
Contro Patient Endogen Pathoge Pathoge	<u>ols</u> Extraction Contro ous Positive Contr n Positive Control n Negative Control	l <sup>1</sup> PASS ol <sup>1</sup> PASS <sup>2</sup> PASS <sup>3</sup> PASS	(1) Endogenous control confirms sample collection, DNA/RNA extraction, and assay (2) Positive control is synthetic inactive pathogen (3) Negative Control contains primers, probe, and enzymes with no DNA/RNA templat (4) A "Detected" result indicates the presence of a pathogen (99.99% confidence) cutoff.			on, and assay enzyme activity NA/RNA template % confidence) above the assay	
Test Pe	rformed	Lab Result (Qualitative Low/M	: (4) edium/High)	DNA Copy Number	/	Comme	ents
SUMMAR	RY Wound/Dermatology	/ Panel			Co	ollection Type: Not Specifi	ed
Enteroc	occus spp.	DETECTED	- LOW	1.00E+02			
Escheri	chia coli	DETECTED	- HIGH	1.00E+04	[08	/23/22] Shares homology with Shig	ella spp.
Staphyl	ococcus epidermidis	DETECTED - N	MEDIUM	1.00E+03			

Class (Gene Name)	Lab Result (Qualitative)	Resistance Gene Targets Identified	Associated Resistances (Antibiotics to Avoid)	
SLIMMARY Wound Antibiotic Resistance			Collection Type: Not Specified	Ī

Markers			
mecA (Methicillin/Oxacillin resistance)	DETECTED	mecA	[08/23/22]Oxacillin

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

FOR MORE WOULD/INFECTION PANEL SAMPLE LAB REPORTS, CLICK HERE

**streamline** 

Page: 1



Sex: M

2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Laboratory

Suite 207

Name:

Address:

**PATIENT INFORMATION** Streamline Pathology Name: PATIENT, TEST DOB: 6/30/1980 2868 Acton Road (205)332-3160 Phone: BIRMINGHAM, AL 35243 Address: 2868 ACTON ROAD Provider: Test, Doctor MD VESTAVIA, AL35243

#### **SPECIMEN INFORMATION**

Lab Accession Number:	S-08232210336
Date Collected:	08/22/22 12:00:00 AM
Date Accessioned:	08/23/2022
Date Reported:	08/23/2022
Faxed to:	18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Wound/Dermatology Panel			Collection Type: Not Specified
Acinetobacter baumannii	Not Detected	Not Detected	
Bacteroides fragilis	Not Detected	Not Detected	
Citrobacter freundii/braakii	Not Detected	Not Detected	
Citrobacter koseri	Not Detected	Not Detected	
Enterobacter cloacae	Not Detected	Not Detected	
Enterococcus spp.	DETECTED - LOW	1.00E+02	
Escherichia coli	DETECTED - HIGH	1.00E+04	[08/23/22] Shares homology with Shigella spp.
Klebsiella aerogenes	Not Detected	Not Detected	
Klebsiella oxytoca/michiganesis	Not Detected	Not Detected	
Klebsiella pneumoniae	Not Detected	Not Detected	
Morganella morganii	Not Detected	Not Detected	
Proteus mirabilis	Not Detected	Not Detected	
Pseudomonas aeruginosa	Not Detected	Not Detected	
Staphylococcus aureus	Not Detected	Not Detected	
Staphylococcus epidermidis	DETECTED - MEDIUM	1.00E+03	
Staphylococcus saprophyticus	Not Detected	Not Detected	
Streptococcus pyogenes (Group A)	Not Detected	Not Detected	
Varicella Zoster virus	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE WOULD/INFECTION PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com



Page: 2



Class (Gene Name)

#### **Molecular PCR Summary Lab Report**

2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

**Associated Resistances** 

CLINIC INFORMATION		PATIENT INFORMATION				SPECIMEN INFORMATION	
Name:	Streamline Pathology Laboratory	Name:	PATIENT, TEST			Lab Accession Number:	S-08232210336
Address:	2868 Acton Road Suite 207	DOB: Phone:	6/30/1980 (205)332-3160	Sex: I	M	Date Collected: Date Accessioned:	08/22/22 12:00:00 AM 08/23/2022
BIRMINGHAM, AL 35243 Provider: Test, Doctor MD		Address:	2868 ACTON ROAD VESTAVIA, AL35243			Date Reported: Faxed to:	08/23/2022 18777966185

**Resistance Gene** 

Lab Result

	(Qualitative)	Targets Identified	(Antibiotics to Avoid)
wound Antibiotic Resistance Ma	arkers		Collection Type: Not Specified
Class A Beta-lactamase (blaKPC)	Not Detected	KPC-2-8,10,11,13-22,24-33	[08/23/22] Carbapenems, Cephalosporins, Penicillins, Beta-lactamase inhibitors, Aztreonam
Class A Beta-lactamase (CTX-M- Group 1)	Not Detected	blaCTX-M-1,3,10,12,15,22,23,28; blaFEC-1	[08/23/22] Cephalosporins, Penicillins, Aztreonam
Class B metallo Beta-lactamase (blaNDM)	Not Detected	NDM (1-21)	[08/23/22] Carbapenems, Cephalosporins, Penicillins, Beta-lactamase inhibitors
Fluoroquinolones	Not Detected	qnrS 1-5,7-9; qnrB Group 1; qnrB Group 5	[08/23/22] Ciprofloxacin, Gemifloxacin, Levofloxacin, Moxifloxacin, Norfloxacin, Ofloxacin
mecA (Methicillin/Oxacillin resistance)	DETECTED	mecA	[08/23/22] Oxacillin
Sulfonamides	Not Detected	sul1; sul2; sul3	[08/23/22] Sulfadiazine, Sulfamethizole, Sulfamethoxazole, Sulfasalazine, Sulfisoxazole
Trimethoprim	Not Detected	dfrA1; dfrA5; dfrA11; dfrA17	[08/23/22] Primsol
vanA (Vancomycin)	Not Detected	vanA	[08/23/22] Vancomycin
vanB (Vancomycin)	Not Detected	vanB	[08/23/22] Vancomycin

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 3

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive.

In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

FOR MORE WOULD/INFECTION PANEL SAMPLE LAB REPORTS, CLICK HERE



## STI PANEL

## PATHOGENS

- Atopobium vaginae
- Chlamydia trachomatis
- Gardnerella vaginalis
- Haemophilus ducreyi
- HHV-1 (Herpes Simplex)
- HHV-2 (Herpes Simplex)
- Neisseria gonorrhoeae
- Treponema pallidum
- Trichonomas vaginalis

For panel offerings, please visit streamlinesci.com



## STI SAMPLE COLLECTION

If submitting a urine sample: Collect a urine sample from the patient in urine collection cup with lid-integrated transfer device: Replace and tighten the collection cup lid, and transfer the urine sample into a vacutainer tube with preservative:

- 1. Remove the sticker from the lid to access the integrated transfer device.
- 2. Insert the vacutainer tube vertically into the transfer device, puncturing the rubber septum and allowing the tube to fill with urine.
- 3. Remove the vacutainer tube and discard the collection cup. Mix urine and preservative by gently inverting the tube 8-10 times.

If submitting a swab sample: Collect a swab sample from the patient with a synthetic flocked swab by following the instructions below depending on the specimen source:



#### Vaginal swab:

- 1. Remove the swab applicator and collect a specimen by rotating the swab against the wall of the vaginal canal several times for 20-30 seconds.
- 2. Withdraw the swab without touching the vaginal surface.

#### Urethral swab:

- 1. Gently insert the swab into the urethra (1-2 cm for women, 2-4 cm for men).
- 2. Rotate the swab in one direction for a minimum of 10 seconds.
- 3. Withdraw the swab.

\*If an ulcer is the intended specimen source, please refer to wound specimen collection instructions

#### Valid specimen type for Chlamydia trachomatis/Neisseria gonorrhoeae/ Trichomonas vaginalis testing only:

#### Throat swab (Oropharyngeal swab):

- 1. Insert the swab into the posterior pharynx and tonsillar areas.
- 2. Rub the swab over both tonsillar pillars and posterior oropharyngeal avoid touching the tongue, teeth, and gums.

#### Rectal swab:

- 1. Insert a sterile swab approximately 2.5 cm into the anal canal.
- 2. Move the swab from side to side in the anal canal. Allow the swab to remain 10-30 seconds for the absorption of organisms onto the swab.
- 3. Remove the swab and insert it into a vial containing 1-3ml of transport media.

#### Disclaimer/Collection - Note:

PCR tests can be used to confirm a suspicion of the presence of sexually transmitted infections in the routine clinical setting. However, in cases of rape or sexual abuse in children under the age of 15, confirmatory testing (along with the appropriate chain of custody as outlined by the CDC) should be performed as recommended by the CDC.

Neisseria gonorrhoeae assay may also detect Neisseria meningitidis due to sequence homology; therefore, a positive result may be interpreted as either Neisseria gonorrhoeae or Neisseria meningitidis for oral STI screening. The patient should be treated based on the clinical presentation of symptoms.

HHV1 assay may also detect HHV2 due to sequence homology; therefore, a positive result may be interpreted as HHV1 or HHV2 infection. The patient should be treated based on the clinical presentation of symptoms.

When collecting a urethral swab or urine specimen from a male or female patient, the patient should not have urinated for **at least an hour** before the specimen is collected. The first-morning urine specimen is preferred due to a large number of cells usually present.

**NOTE:** Certain organisms are intracellular; therefore, there must be enough human cells present to detect the organism.





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

CLINIC I	NFORMATION	PATIENT	INFORMATION		SPECIMEN INFORMAT	ΓΙΟΝ
Name: Address: Provider:	Streamline Pathology Laboratory 2868 Acton Road Suite 207 BIRMINGHAM, AL 35243 Test, Doctor MD	Name: DOB: Phone: Address:	PATIENT, TEST 6/30/1980 (205)332-3160 2868 ACTON ROAD VESTAVIA, AL35243	Sex: M	Lab Accession Number: Date Collected: Date Accessioned: Date Reported: Faxed to:	S-08242210101 08/22/22 12:00:00 AM 08/24/2022 08/24/2022 18777966185
Contr Patient Endogen Pathoge Pathoge	<u>ols</u> Extraction Control <sup>1</sup> ious Positive Control <sup>2</sup> in Positive Control <sup>2</sup> in Negative Control <sup>3</sup>	PASS PASS PASS PASS	<ol> <li>Endogenous control a</li> <li>Positive control is</li> <li>Negative Control con</li> <li>A "Detected" result cutoff.</li> </ol>	confirms sampl synthetic ina ntains primers indicates the	e collection, DNA/RNA extracti active pathogen 3, probe, and enzymes with no D 2 presence of a pathogen (99.99	on, and assay enzyme activity NA/RNA template % confidence) above the assay

PASS cutoff.

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments

SUMMARY Sexually Transmittee	I Infection Pathogens	Collection Type: Swab	
Atopobium vaginae	DETECTED - HIGH	1.00E+06	
HHV-1 (Herpes Simplex Virus)	DETECTED - MEDIUM	1.00E+03	[08/24/22] HHV-1 assay may also detect HHV-2 due to sequence homology; therefore, a positive result may be interpreted as HHV-1 or HHV-2 infection. The patient should be treated based on the clinical presentation of symptoms.

Clinic Name: Streamline Pathology Laboratory

Lab Director: Ty Thomas, MD

Page: 1

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

FOR MORE STI PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

CLINIC INFORMATION		PATIENT INFORMATION			SPECIMEN INFORMATION	
Name:	Streamline Pathology Laboratory	Name:	PATIENT, TEST	_	Lab Accession Number:	S-08242210101
Address:	2868 Acton Road Suite 207	DOB: Phone:	6/30/1980 (205)332-3160	Sex: M	Date Collected: Date Accessioned:	08/22/22 12:00:00 AM 08/24/2022
BIRMINGHAM, AL 35243 Provider: Test, Doctor MD		Address:	dress: 2868 ACTON ROAD VESTAVIA, AL35243		Date Reported: Faxed to:	08/24/2022 18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Sexually Transmitted Infection	Pathogens		Collection Type: Swab
Atopobium vaginae	DETECTED - HIGH	1.00E+06	
Chlamydia trachomatis	Not Detected	Not Detected	
Gardnerella vaginalis	Not Detected	Not Detected	
Haemophilus ducreyi	Not Detected	Not Detected	
HHV-1 (Herpes Simplex Virus)	DETECTED - MEDIUM	1.00E+03	[08/24/22] HHV-1 assay may also detect HHV-2 due to sequence homology; therefore, a positive result may be interpreted as HHV-1 or HHV-2 infection. The patient should be treated based on the clinical presentation of symptoms.
HHV-2 (Herpes Simplex Virus)	Not Detected	Not Detected	
Neisseria gonorrhoeae	Not Detected	Not Detected	
Treponema pallidum	Not Detected	Not Detected	
Trichomonas vaginalis	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 2

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE STI PANEL SAMPLE LAB REPORTS, CLICK HERE



## VAGINITIS PANEL

## PATHOGENS

- Atopobium vaginae
- Bacteroides fragilis
- BVAB-2
- Candida albicans
- Candida dubliniensis
- Candida glabrata
- Candida krusei
- Candida lusitaniae
- Candida parapsilosis
- Candida tropicalis
- Chlamydia trachomatis
- Enterococcus spp.
- Escherichia coli
- Gardnerella vaginalis
- Haemophilus ducreyi
- HHV-1 (Herpes Simplex)
- HHV-2 (Herpes Simplex)

- Lactobacillus crispatus
- Lactobacillus gasseri
- Lactobacillus iners
- Lactobacillus jensenii
- Megasphaera Type 1
- Megasphaera Type 2
- Mobiluncus curtisii
- Mobiluncus mulieris
- Mycoplasma genitalium
- Mycoplasma hominis
- Neisseria gonorrhoeae
- Prevotella bivia
- Staphylococcus aureus
- Streptococcus agalactiae (Group B)
- Treponema pallidum
- Trichonomas vaginalis
- Ureaplasma urealyticum



## VAGINITIS SAMPLE COLLECTION



Collect a vaginal sample from the patient using a synthetic flocked swab by inserting the swab into the vagina and turn the swab 3 times.

#### **Disclaimer/Collection Note:**

PCR tests can be used to confirm a suspicion of the presence of sexually transmitted infections in the routine clinical setting. However, in cases of rape or sexual abuse in children under the age of 15, confirmatory testing (along with the appropriate chain of custody as outlined by the CDC) should be performed as recommended by the CDC.

HHV1 assay may also detect HHV2 due to sequence homology; therefore, a positive result may be interpreted as HHV1 or HHV2 infection. The patient should be treated based on the clinical presentation of symptoms.

**NOTE:** Certain organisms are intracellular; therefore, there must be enough human cells present to detect the organism.





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

	NFORMATION	PATIENT	INFORMATION			SPECIMEN INFORMAT	ION
Name: Address:	Streamline Pathology Laboratory 2868 Acton Road Suite 207	Name: DOB: Phone:	PATIENT, TEST 6/30/1980 (205)332-3160	Sex:	М	Lab Accession Number: Date Collected: Date Accessioned:	S-08232210335 08/22/22 12:00:00 AM 08/23/2022
Provider:	Test, Doctor MD	Address:	2868 ACTON ROAD VESTAVIA, AL35243			Date Reported: Faxed to:	08/23/2022 18777966185
			(1) Endogonous control (	onfirm		o collection DNA/BNA extracti	on and accay ongume activity

<u>Controls</u>			
Patient Extraction Control <sup>1</sup>	PASS	(2)	
Endogenous Positive Control <sup>1</sup>	PASS	(3)	
Pathogen Positive Control <sup>2</sup>	PASS	(4)	
Pathogen Negative Control 3	PASS	cut	

Positive control is synthetic inactive pathogen

Negative Control contains primers, probe, and enzymes with no DNA/RNA template

A "Detected" result indicates the presence of a pathogen (99.99% confidence) above the assay off.

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
SUMMARY Vaginitis Pathoge	ens		Collection Type: Not Specified
Candida dubliniensis	DETECTED - HIGH	1.00E+05	
Candida tropicalis	DETECTED - LOW	1.00E+02	
Lactobacillus crispatus	DETECTED - MEDIUM	8.90E+02	
Lactobacillus gasseri	DETECTED - HIGH	8.37E+03	
Lactobacillus iners	DETECTED - LOW	1.83E+02	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE VAGINITIS PANEL SAMPLE LAB REPORTS, CLICK HERE

PCR CLIENT SERVICES MANUAL • www.streamlinesci.com



Page: 1



Μ

2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Name:	Streamline Pathology Laboratory
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243
Provider:	Test, Doctor MD

PATIENT INFORMATION				
Name:	PATIENT, TEST			
DOB:	6/30/1980	Sex:		
Phone:	(205)332-3160			
Address:	2868 ACTON ROAD VESTAVIA, AL35243			

SPECIMEN INFORMATION				
Lab Accession Number:	S-08232210335			
Date Collected:	08/22/22 12:00:00 AM			
Date Accessioned:	08/23/2022			
Date Reported:	08/23/2022			
Faxed to:	18777966185			

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Vaginitis Pathogens			Collection Type: Not Specified
Atopobium vaginae	Not Detected	Not Detected	
Bacteroides fragilis	Not Detected	Not Detected	
BVAB-2	Not Detected	Not Detected	
Candida albicans	Not Detected	Not Detected	
Candida dubliniensis	DETECTED - HIGH	1.00E+05	
Candida glabrata	Not Detected	Not Detected	
Candida krusei	Not Detected	Not Detected	
Candida lusitaniae	Not Detected	Not Detected	
Candida parapsilosis	Not Detected	Not Detected	
Candida tropicalis	DETECTED - LOW	1.00E+02	
Chlamydia trachomatis	Not Detected	Not Detected	
Enterococcus spp.	Not Detected	Not Detected	
Escherichia coli	Not Detected	Not Detected	
Gardnerella vaginalis	Not Detected	Not Detected	
Haemophilus ducreyi	Not Detected	Not Detected	
HHV-1 (Herpes Simplex Virus)	Not Detected	Not Detected	
HHV-2 (Herpes Simplex Virus)	Not Detected	Not Detected	
Lactobacillus crispatus	DETECTED - MEDIUM	8.90E+02	
Lactobacillus gasseri	DETECTED - HIGH	8.37E+03	
Lactobacillus iners	DETECTED - LOW	1.83E+02	
Lactobacillus jensenii	Not Detected	Not Detected	
Megasphaera Type 1	Not Detected	Not Detected	
Megasphaera Type 2	Not Detected	Not Detected	
Mobiluncus curtisii	Not Detected	Not Detected	
Mobiluncus mulieris	Not Detected	Not Detected	
Mycoplasma genitalium	Not Detected	Not Detected	
Mycoplasma hominis	Not Detected	Not Detected	
Neisseria gonorrhoeae	Not Detected	Not Detected	
Prevotella bivia	Not Detected	Not Detected	
Staphylococcus aureus	Not Detected	Not Detected	
Streptococcus agalactiae (GBS)	Not Detected	Not Detected	
Treponema pallidum	Not Detected	Not Detected	
Trichomonas vaginalis	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE VAGINITIS PANEL SAMPLE LAB REPORTS, CLICK HERE

. . . . . . . . . . . . . . . . .



Page: 2



2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

CLINIC IN	IFORMATION	PATIENT	INFORMATION		SPECIMEN INFORMAT	ION
Name:	Streamline Pathology Laboratory	Name:	PATIENT, TEST		Lab Accession Number:	S-08232210335
Address.	2868 Acton Boad	DOB:	6/30/1980	Sex: M	Date Collected:	08/22/22 12:00:00 AM
Address.	Suite 207	Phone:	(205)332-3160		Date Accessioned:	08/23/2022
<b>_</b>	BIRMINGHAM, AL 35243	Address:	2868 ACTON ROAD		Date Reported:	08/23/2022
Provider:	Test, Doctor MD		VESTAVIA, AL35243		Faxed to:	18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Ureaplasma urealyticum	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 3

streamline

scientific

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive every sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detectd, resulting in a false negative result.

FOR MORE VAGINITIS PANEL SAMPLE LAB REPORTS, CLICK HERE



## GASTROINTESTINAL PANEL

## PATHOGENS

- Adenovirus
- Astrovirus
- Campylobacter (coli, jejuni, lari)
- Clostidium difficile
  - c. difficile toxin A
  - c. difficile toxin B
- Escherichia coli (VTEC)
  - Shiga-toxin 1
  - Shiga-toxin 2
  - E.Coli )157
- Norovirus GI
- Norovirus GII
- Rotavirus
- Salmonella spp.
- Sapovirus
- Shigella spp.
- Yersinia enterocolitica

For panel offerings, please visit streamlinesci.com



## GASTROINTESTINAL SAMPLE COLLECTION



Collect a gastrointestinal stool specimen sample or a rectal swab from the patient using a swab with Cary Blair buffer (ex: Copan Eswab). Use the Cary Blair buffer specimen collection swab for the transport of either specimen swab or rectal swab by using the following instructions:

- 1. Remove the swab and transport tube from the packaging. Do not contaminate.
- 2. For stool specimen swabs, carefully collect a portion of the specimen by placing the swab tip directly into the specimen.
- 3. For rectal swabs, insert the tip of the swab approximately 1 inch beyond the anal sphincter.
- 4. Carefully rotate the swab to sample the anal crypts, then withdraw the swab.





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

	NFORMATION	PATIENT	INFORMATION			SPECIMEN INFORMAT	ION
Name: Address: Provider:	Streamline Pathology Laboratory 2868 Acton Road Suite 207 BIRMINGHAM, AL 35243 Test, Doctor MD	Name: DOB: Phone: Address:	PATIENT, TEST 6/30/1980 (205)332-3160 2868 ACTON ROAD VESTAVIA, AL35243	Sex: N	И	Lab Accession Number: Date Collected: Date Accessioned: Date Reported: Faxed to:	S-08232210337 08/22/22 12:00:00 AM 08/23/2022 08/24/2022 18777966185
Contro Patient Endogen Pathoge Pathoge	Dls Extraction Control <sup>1</sup> ous Positive Control <sup>1</sup> n Positive Control <sup>2</sup> n Negative Control <sup>3</sup>	PASS PASS PASS PASS	<ol> <li>(1) Endogenous control (2) Positive control is</li> <li>(3) Negative Control co</li> <li>(4) A "Detected" result cutoff.</li> </ol>	confirms s synthetic ntains pri indicates	sampl c ina imers s the	e collection, DNA/RNA extraction ctive pathogen , probe, and enzymes with no Di presence of a pathogen (99.99)	on, and assay enzyme activity NA/RNA template % confidence) above the assay

Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Panel - Bacterial Pathogens		Collection Type: Not Specified, Stool - Fecal
DETECTED - LOW	DETECTED - LOW	
DETECTED - HIGH	1.00E+04	[08/23/22] Salmonella was detected by PCR and is a Notifiable Pathogen per ADPH. Result to be confirmed by lab via bacterial culture before reporting to ADPH. The required time for the ordering provider to notify ADPH is within 5 days from receipt of report.
SUMMARY Gastrointestinal Panel - Viral Pathogens		Collection Type: Not Specified, Stool - Fecal
	Lab Result (4) (Qualitative Low/Medium/High) Panel - Bacterial Pathogens DETECTED - LOW DETECTED - HIGH Panel - Viral Pathogens	Lab Result (4) DNA Copy (Qualitative Low/Medium/High) Number Panel - Bacterial Pathogens DETECTED - LOW DETECTED - LOW DETECTED - HIGH 1.00E+04 Panel - Viral Pathogens

Clinic Name: Streamline Pathology Laboratory

Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 1

streamline

scientific

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive. In addition, while these assays are very sensitive may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false positive.

FOR MORE GASTROINTESTINAL PANEL SAMPLE LAB REPORTS, CLICK HERE





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### **CLINIC INFORMATION**

Name:	Streamline Pathology Laboratory
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243
Provider:	Test, Doctor MD

PATIENT INFORMATION					
Name:	PATIENT, TEST				
DOB:	6/30/1980	Sex:	М		
Phone:	(205)332-3160				
Address:	2868 ACTON ROAD VESTAVIA, AL35243				

### SPECIMEN INFORMATION

Lab Accession Number:	S-08232210337
Date Collected:	08/22/22 12:00:00 AM
Date Accessioned:	08/23/2022
Date Reported:	08/24/2022
Faxed to:	18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Gastrointestinal Panel - Bacte	erial Pathogens		Collection Type: Not Specified, Stool - Fecal
Campylobacter coli/jejuni/lari	Not Detected	Not Detected	
Clostridium difficile	Not Detected	Not Detected	
Toxin A (tcdA)	Not Detected	Not Detected	-
Toxin B (tcdB)	Not Detected	Not Detected	
Escherichia coli (VTEC)	DETECTED - LOW	DETECTED - LOW	
O157	Not Detected	Not Detected	
STX1	Not Detected	Not Detected	
STX2	Not Detected	Not Detected	
Salmonella spp.	DETECTED - HIGH	1.00E+04	[08/23/22] Salmonella was detected by PCR and is a Notifiable Pathogen per ADPH. Result to be confirmed by lab via bacterial culture before reporting to ADPH. The required time for the ordering provider to notify ADPH is within 5 days from receipt of report.
Shigella spp.	Not Detected	Not Detected	
Yersinia enterocolitica	Not Detected	Not Detected	

Gastrointestinal Panel - Viral Pathogens			Collection Type: Not Specified, Stool - Fecal
Adenovirus	Not Detected	Not Detected	
Astrovirus	Not Detected	Not Detected	
Norovirus GI	Not Detected	Not Detected	
Norovirus GII	Not Detected	Not Detected	
Rotavirus	Not Detected	Not Detected	
Sapovirus	Not Detected	Not Detected	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 2

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE GASTROINTESTINAL PANEL SAMPLE LAB REPORTS, CLICK HERE



## FUNGAL INFECTION PANEL

## PATHOGENS

- Alternaria spp.
- Aspergillus spp.
- Fusarium spp.
- Scytalidium dimidiatum
- Sarocladium strictum
- Candida albicans
- Candida glabrata
- Candida krusei
- Candida parapsilosis

- Candida tropicalis
- Cryptococcus spp.
- Malassezia spp.
- Meyerozyma guillermondii
- Trichophyton anthropophilic spp.

. . . . . . . . . . . . . . . .

. . . . . . . . . . . . . . . . .

- Trichophyton zoophilic spp.
- Mircosporum canis

## **BACTERIA ADD ON**

• Pseudomonas aeruginosa

## ABX RESISTANCE MARKER

. . . . .

• Methicillin/Oxacillin (mecA)



## FUNGAL SPECIMEN COLLECTION

## ACCEPTABLE SPECIMENS:

Nail clippings and skin scrapings in a dry collection tube. If a wound is infected, the area should be collected by a swab in liquid amies media.

### DIRECTIONS FOR COLLECTING A NAIL CLIPPING

- 1. Wipe the nail collection site with 70% isopropyl alcohol.
- 2. Debride and discard nail clippings.
- 3. Obtain specimen from the most proximal area of nail and hyponychium. Minimum specimen amount size of nail and subungual debris is 3mm to 6mm (small pieces to obtain this size are preferred).
- 4. Use a curette to obtain any additional subungual debris.
- 5. Place the dry nail sample and debris into a dry sterile tube.

### DIRECTIONS FOR A SKIN SCRAPING

- 1. Remove any traces of skin products, medications, or surface contaminants by wiping the area with a 70% isopropyl alcohol wipe.
- 2. Choose the best area to scrape by determining where fungal growth is most active.
- 3. Scrape the skin using a scalpel held at a blunt angle into dry collection tube. The greater amount of specimen, the better the result.

### DIRECTIONS FOR A WOUND COLLECTION

- 1. Cleanse and debride the infected area with 70% isopropyl alcohol using sterile gauze. Saline can be used if there is an open wound.
- 2. Exudate and brushing of the base of the infected area (including advancing margins) should be collected using the swab. Without contaminating the swab, place the swab in the media.



2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

the assay

CLINIC I	NFORMATION	PATIENT	INFORMATION			SPECIMEN INFORMAT	ΓΙΟΝ
Name: Address: Provider:	Streamline Pathology Laboratory 2868 Acton Road Suite 207 BIRMINGHAM, AL 35243 Test Doctor MD	Name: DOB: Phone: Address:	PATIENT, TEST 6/30/1980 (205)332-3160 2868 ACTON ROAD	Sex:	М	Lab Accession Number: Date Collected: Date Accessioned: Date Reported:	S-08232210339 08/22/22 12:00:00 AM 08/23/2022 08/23/2022
						Faxed to:	18777966185
Contro Patient Endogen	<b>ols</b> Extraction Control <sup>1</sup> Lous Positive Control <sup>1</sup>	PASS PASS	<ol> <li>Endogenous control of</li> <li>Positive control is</li> <li>Negative Control control</li> </ol>	synthet	sampl ic ina rimers	e collection, DNA/RNA extracti ctive pathogen , probe, and enzymes with no D	on, and assay enzyme activity NA/RNA template
Pathoge	n Positive Control <sup>2</sup>	PASS	(4) A "Detected" result	indicat	og tho	presence of a pathogen (99 99	<pre>% confidence) above the assau</pre>

SUMMARY Fungal Infection Panel			Collection Type: Swab			
Test Performed	Lab Result (Qualitative Low/M	: (4) edium/High)	DNA Copy Number		Comments	
Pathogen Positive Control Pathogen Negative Control	<sup>2</sup> PASS <sup>3</sup> PASS	(4) A "Detect cutoff.	ted" result indicate	s the presence o	of a pathogen (99.99% confidence) al	DOVE

Sommart rungar intection ranei			Collection Type: Swab
Meyerozyma guilliermondii	DETECTED - MEDIUM	1.00E+03	
Trichosporon spp.	DETECTED - HIGH	1.00E+06	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

#### Lab Director: Ty Thomas, MD

Page: 1

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE FUNGAL INFECTION PANEL SAMPLE LAB REPORTS, CLICK HERE





2868 Acton Road, Suite 207 Birmingham, AL, 35243 Phone: (855) 319-4459 CLIA ID: 01-D2074949

#### CLINIC INFORMATION

Name:	Streamline Pathology Laboratory
Address:	2868 Acton Road Suite 207 BIRMINGHAM, AL 35243
Provider:	Test, Doctor MD

PATIENT INFORMATION				
Name:	PATIENT, TEST			
DOB:	6/30/1980	Sex:	Μ	
Phone:	(205)332-3160			
Address:	2868 ACTON ROAD VESTAVIA, AL35243			

# SPECIMEN INFORMATIONLab Accession Number:S-08232210339Date Collected:08/22/22 12:00:00 AMDate Accessioned:08/23/2022Date Reported:08/23/2022Faxed to:18777966185

Test Performed	Lab Result (4) (Qualitative Low/Medium/High)	DNA Copy Number	Comments
Fungal Infection Panel			Collection Type: Swab
Alternaria spp.	Not Detected	Not Detected	
Aspergillus spp.	Not Detected	Not Detected	
Candida albicans	Not Detected	Not Detected	
Candida glabrata	Not Detected	Not Detected	
Candida krusei	Not Detected	Not Detected	
Candida parapsilosis	Not Detected	Not Detected	
Candida tropicalis	Not Detected	Not Detected	
Cryptococcus spp.	Not Detected	Not Detected	
Curvularia spp.	Not Detected	Not Detected	
Epidermophyton floccosum	Not Detected	Not Detected	
Fusarium spp.	Not Detected	Not Detected	
Malassezia spp.	Not Detected	Not Detected	
Meyerozyma guilliermondii	DETECTED - MEDIUM	1.00E+03	
Microsporum canis	Not Detected	Not Detected	
Sarocladium strictum	Not Detected	Not Detected	
Scytalidium dimidiatum	Not Detected	Not Detected	
Trichophyton anthropophilic spp.	Not Detected	Not Detected	[08/23/22] Targets: T. rubrum, T. violaceum, T. soudanense, T. concentricum and other species
Trichophyton zoophilic spp.	Not Detected	Not Detected	[08/23/22] Targets: T. mentagrophytes, T. interdigitale (anthropophilic but related to T. mentagrophytes), T. tonsurans (anthropophilic but related to T. equinum), T. verrucosum, T. simii, T. benhamiae, T. equinum, and other species
Trichosporon spp.	DETECTED - HIGH	1.00E+06	

Clinic Name: Streamline Pathology Laboratory

#### Patient name: PATIENT, TEST

Lab Director: Ty Thomas, MD

Page: 2

Test results must be evaluated with clinical symptoms to diagnose disease. All tests established and validated by Laboratory and not FDA approved unless otherwise indicated. Methodology statement: Real-time PCR assays are designed to detect pathogens with clinical significance and analytical sensitivity and specificity greater than 99%. Limitations: While these assays are very sensitive and specific, theoretically these assays could detect pathogens not listed, resulting in a false positive. In addition, while these assays are very specific, there may be target pathogen sequences with unknown sequence variability which may not be detected, resulting in a false negative result.

#### FOR MORE FUNGAL INFECTION PANEL SAMPLE LAB REPORTS, CLICK HERE



